



NEANN FAB



First Aid Backpack Layout Guide

PLEASE READ OR WARRANTY MAY BE VOID



The following booklet is a suggested layout guide for setting up the NEANN FAB. It should be used as a guideline only.





MAIN COMPARTMENT





FRONT COMPARTMENT



Velcro Attachments
for
IV And ETT Kits



Pneumocath

IV Giving
Sets

IV Fluids

Intraosseous

Triangular Bandages



IV ACCESS KIT





ENDOTRACHEAL KIT





CLEANING INSTRUCTIONS

All fabrics used in the FAB are designed to comply with the **Australian Standards on Laundry Practice (AS 4146-1994)** for the removal and killing of HIV, Hepatitis B and Vegetative Organisms.

A cleaning program that complies with the Australian Standards is available directly from RAPP Australia Pty Ltd for cleaning blood, vomit and other substances off the FAB. For those wishing to undertake in-house cleaning (at you own risk, including risk to warranty), some basic guidelines are listed below:

- *The cleaner in keeping with normal practices should put on personal protective equipment such as gloves, safety glasses, face masks, etc..*
- *Small areas of contamination:*
 - Blood - Soak by placing a wet sponge on stain immediately after contact. If available, an antibacterial solution (Virkon) should be added to the water. Repeat as necessary, then gently clean off residual stains with light spray of Nifty. Allow to dry thoroughly before use.*
 - Vomit - Gently sponge with hot water immediately after contact. If available, an antibacterial solution (Virkon) should be added to the water. Repeat as necessary, then gently clean off residual stains with light spray of Nifty. Allow to dry thoroughly before use.*
- *Larger areas of contamination:*
 - Initially - Remove resuscitator components carefully . Remove stiffening corflute from effected panels.*
 - Blood - Soak by placing effected panels in cold water. An antibacterial solution (Virkon) should be added to the water. Let effected part soak for 20 minutes. Using Nifty, gently sponge of residual stains. Allow to dry thoroughly before use .*
 - Vomit - Soak by placing in hot water. An antibacterial solution (Virkon) should be added to the water. Using Nifty, gently sponge of residual stains. Allow to dry thoroughly before use.*

In very heavy saturation, the FAB can be dismantled, soaked in cold water with an antibacterial solution (Virkon) for 2 hours, rinsed and then emersed again and brought up to a temperature of 90 degrees Celsius and maintained for 10 minutes. Using Nifty, gently sponge of residual stains. Dry thoroughly before reassembly and use.

Failure to comply with proper cleaning instructions may damage material and void warranty



WARRANTY

With 25 years of manufacturing experience and a commitment to using only the highest quality materials, Neann's range of medical bags are guaranteed against faulty workmanship for the life of the bag. This means if the bag fails due to a manufacturing fault, it will be repaired or replaced at the manufacturers option. This of course excludes fair wear and tear, or misuse of the bag against manufacturers guidelines.

NEANN

EMERGENCY CARE PRODUCTS

By

RAPP Australia Pty Ltd
P.O. Box
Lara
Victoria 3212

Phone: (03) 5284 0222

Fax: (03) 5282 4892



AUSTRALIAN MADE®



NEANN FAB – FIRST AID BACKPACK **USER SURVEY**

To help RAPP improve the product & training information supplied, please complete the following survey by circling the appropriate number next to the question. Additional comments are most welcome.

Name:

Station:

Contact Phone number:

Very poor
Poor
Average
Good
Very Good
Excellent

1. Was the powerpoint training package adequate? 1 2 3 4 5

Comment:

2. Was this information booklet adequate? 1 2 3 4 5

Comment:

3. Are you pleased with the manufacturing quality? 1 2 3 4 5

Comment:

4. Does all the equipment fit properly into the kit? 1 2 3 4 5

Comment:

5. Is equipment easy to access? 1 2 3 4 5

Comment:

6. Does the kit work well in the field? 1 2 3 4 5

Comment:

7. Are you pleased with the kit overall? 1 2 3 4 5

Comment:

Once complete, please fax to the RAPP Office on (03) 5282 4892