

LSP Aspirator Unit

Information

LSP Aspirator Unit Information



This Photo shows the standard venturi, operation button and threaded inlet.

This photo shows a comparison of the standard venturi (Top) and the newer venturi and redirector (Bottom).



The new venturi is ported so as to maintain the standard amount of negative pressure production and also accommodate the longer exhaust/condenser tube.

This photo shows the condenser tube attached. The condenser tube prevents body fluids being ejected from the aspirator in a mist.



The tube allows for 1/. Condensing of the vapour and 2/. management of the exhaust.