

RES-Q-VAC[®]

Features and Benefits

Features

Benefits

Safety

Full Stop Protection[®]

Full Stop Protection
0.22 Micron Filter

Safeguards personnel from pathogens including:

Flu, HIV/AIDS, TB, Hepatitis, and SARS
Meets CDC SARS and OSHA Regulations
Perfect for Isolated Patients
Provides overflow protection

Performance

Maximum Vacuum- >600mm hg

Power to clear vomit, clotting blood, and mucus
Exceeds hospital & ambulance standards

Continuous Suctioning

Suctions until fluids clear tubing
Less pumping required

Manual Power

Assures suction when needed
No battery maintenance, no O₂ consumption
A necessity when dealing with bio-chemical warfare

Molded ABS Plastic

Light and Rugged
Operates from 0-140 degrees F
Gama sterilized

Easy tip control

Ease of Use

Hand Held

No long tubing to handle or maintain

Mechanically Operated

No need for electric or batteries
Works when power fails, or isn't available

Easy to Pull Lever

Will not fatigue user
Minimal effort required
Achieve high sustained air flow

Angled Tubing Yankauer

Allows clear view of airway

Canister Design

Attaches simply to pump
Easy to monitor filling
Can be capped for sampling or disposal

Reusable Pump

Exceeds ISO manual pump standards

RES-Q-VAC®

Features and Benefits

Features

Benefits

Convenience

Size

Pump 4.5" x 7.1" x 2.1"

Easy to use in restrictive situations
Occupies minimal storage space

Weight

Pump 5 oz, Canister 2 oz.

Easy to hold and control

Full Stop Protection® Filter 1.5 oz.

Complies with CDC and OSHA guidelines

Choice of Suction Tips

Soft wide bore

Rigid yankauer

E.T. Adapter

Pediatric Catheters

Adapts to various clinical situations
First Responder through MD requirements

Canister Paks, Sterile

No cross-contamination
No lost parts
Sampling cap and label provided

RMS Medical Products

Proudly offers quality medical products:

DENTAL-EVAC®

FREEDOM60® Infusion System

Gyneco Thermal Cautery Unit

Masterson Endometrial Biopsy Kit

RMS Medical Products

24 Carpenter Road

Chester, NY 10918

Phone 845.469.2042

Fax 845.469.5518

www.rmsmedicalproducts.com