

# ONLY ONE







# air-Q<sup>®</sup>sp



## The only self-pressurising mask cuff

The **air-Qsp** from Mercury Medical provides many advantages in life-saving emergencies, including:

-  **Eliminates mask cuff over-inflation**  
The cyclical lowering in intra-cuff pressure may assist in diminishing complications such as Mucosal and nerve trauma
-  **No inflation line necessary**  
With the new patented self-pressurising, disposable air-Qsp, there is no need to inflate the mask; it self-pressurises to the proper inflation pressure automatically
-  **Eliminates the extra step of inflating the mask cuff**  
The new design allows Positive Pressure Ventilation to self-pressurise the mask cuff. The increase in cuff seal pressure occurs at the exact time you need it, during the upstroke of ventilation  
  
On exhalation, the cuff decompresses to the level of the PEEP
-  **Facilitates intubation**  
The SP version can be used as an 'everyday' airway that facilitates intubation

## 6 SIZES AVAILABLE

**SAFE • SIMPLE • SELF-PRESSURISING  
ELIMINATES OVER-INFLATION**



is a division of RAPP Australia Pty Ltd

Pre-eminent supplier of pre-hospital emergency and rescue equipment for First Responder, Ambulance and Fire Services since 1983

Provider of ancillary emergency products for all emergency situations including military, mining, diving and other recreational activities

### Shared Goals

We share with almost all of our customers these three foundation goals:

- Effectiveness
- Reliability and
- Safety of Use

We make it our business to help our customers to achieve their goals

### Shared Achievement

We successfully work with our customers to ensure that they can provide the most appropriate equipment and supplies to their personnel, both professional and volunteer, to assist them to complete their tasks successfully and ensure they return safely

### Future Commitment

To maintain our leadership position we commit to:

- continuously apply innovation to our existing products,
- work with our customers' design teams to develop new products
- constantly search the world for latest best-practice, and introduce those which fit as soon as practicable



### For more information

Visit

[www.rappaustralia.com.au](http://www.rappaustralia.com.au)

Email

[sales@rappaustralia.com.au](mailto:sales@rappaustralia.com.au)

Call

+61 (0) 3 5284 0222

Postal Address

PO Box 120, Lara, Victoria, Australia 3212