



Mercury Medical®

**ONLY ONE**

# Neo-Tee®

Disposable Infant T-Piece Resuscitator



The first and only disposable infant t-piece resuscitator

The **Neo-Tee** from Mercury Medical provides many advantages in life-saving emergencies, including:

- ✓ **Consistent PIP & PEEP Pressure**  
The Neo-Tee is flow-controlled and pressure limited. It offers the ability to measure a more consistent, targeted Peak Inspiratory Pressure (PIP) and Positive End-Expiratory Pressure (PEEP)
- ✓ **Built-in manometer verifies pressures**  
The portable, built-in Manometer 'on the tee' has an accurate and easy-to-read dial and safety pressure relief
- ✓ **Cost efficient**  
With new Neo-Tee there's no capital equipment to purchase. It is completely disposable, so there is no cleaning required, it is space saving, lightweight and easily transported
- ✓ **Controller**  
Allows clinicians to change airway pressure "in-line" closer to the patient

**FLOW CONTROLLED • PRESSURE LIMITED**  
**CHANGE AIRWAY PRESSURE CLOSER TO THE PATIENT**

**WWW.NEANN.COM.AU**



is a division of RAPP Australia Pty Ltd

Pre-eminent supplier of pre-hospital emergency and rescue equipment for First Responder, Ambulance and Fire Services since 1983

Provider of ancillary emergency products for all emergency situations including military, mining, diving and other recreational activities

### Shared Goals

We share with almost all of our customers these three foundation goals:

- Effectiveness
- Reliability and
- Safety of Use

We make it our business to help our customers to achieve their goals

### Shared Achievement

We successfully work with our customers to ensure that they can provide the most appropriate equipment and supplies to their personnel, both professional and volunteer, to assist them to complete their tasks successfully and ensure they return safely

### Future Commitment

To maintain our leadership position we commit to:

- continuously apply innovation to our existing products,
- work with our customers' design teams to develop new products
- constantly search the world for latest best-practice, and introduce those which fit as soon as practicable



### For more information

Visit

[www.rappaustralia.com.au](http://www.rappaustralia.com.au)

Email

[sales@rappaustralia.com.au](mailto:sales@rappaustralia.com.au)

Call

+61 (0) 3 5284 0222

Postal Address

PO Box 120, Lara, Victoria, Australia 3212