

# Bag Valve Mask (BVM) Options Matrix

| Laerdal Silicone Resuscitators |                      |                      |                      |                       |                      |                      |                      |
|--------------------------------|----------------------|----------------------|----------------------|-----------------------|----------------------|----------------------|----------------------|
|                                | Adult                |                      |                      | Paediatric            |                      | Pre-Term             |                      |
|                                | Patients >25kg       |                      |                      | Patient 2.5kg to 25kg |                      | Patients < 2.5kg     |                      |
|                                | Basic w/o Mask       | Standard             | Complete             | Basic w/o Mask        | Complete             | Basic w/o Mask       | Complete             |
| <b>Carton</b>                  | 50-003-87005033-0000 | 50-003-87005233-0000 | 50-003-87005133-0000 | 50-003-86005033-0000  | 50-003-86005233-0000 | 50-003-85005033-0000 | 50-003-85005133-0000 |
| <b>Compact Case</b>            |                      |                      | 50-003-87005333-0000 |                       | 50-003-86005340-0000 |                      | 50-003-85005333-0000 |
| <b>Display Case</b>            |                      |                      | 50-003-87005533-0000 |                       | 50-003-86005540-0000 |                      | 50-003-85005533-0000 |
| <b>Bag 1600ml</b>              | ✓                    | ✓                    | ✓                    |                       |                      |                      |                      |
| <b>500ml</b>                   |                      |                      |                      | ✓                     | ✓                    |                      |                      |
| <b>240ml</b>                   |                      |                      |                      |                       |                      | ✓                    | ✓                    |
| <b>Pt Valve w/o Relief</b>     | ✓                    | ✓                    | ✓                    |                       |                      |                      |                      |
| <b>w/Relief Valve</b>          |                      |                      |                      | ✓                     | ✓                    | ✓                    | ✓                    |
| <b>Mask 0/0</b>                |                      |                      |                      |                       |                      |                      | ✓                    |
| <b>0/1</b>                     |                      |                      |                      |                       | ✓                    |                      | ✓                    |
| <b>2</b>                       |                      |                      |                      |                       | ✓                    |                      |                      |
| <b>3/4 w/Mask Cover</b>        |                      |                      | ✓                    |                       | ✓                    |                      |                      |
| <b>4/5+ w/Mask Cover</b>       |                      | ✓                    | ✓                    |                       |                      |                      |                      |
| <b>Intake/Reservoir Bag</b>    | ✓                    | ✓                    | ✓                    | ✓                     | ✓                    | ✓                    | ✓                    |
| <b>Reservoir Bag 2600ml</b>    | ✓                    | ✓                    | ✓                    |                       |                      |                      |                      |
| <b>Reservoir Bag 600ml</b>     |                      |                      |                      | ✓                     | ✓                    | ✓                    | ✓                    |