



RESPIRATORY & TRAUMA RANGE



NEANN

NEANN is the medical equipment division of RAPP Australia Pty Ltd. The focus of the division is to design, develop and provide first response, emergency and hospital medical equipment and supplies that are effective, reliable and safe to use.

NEANN is a leader among medical equipment suppliers and our products can be found in use around Australia in most ambulances, hospitals, clinics and emergency services teams. RAPP Australia Pty Ltd holds ongoing contracts with all State Governments in Australia as well as leading emergency agencies worldwide.

The division works in collaboration with clients to ensure all equipment performs as intended. NEANN is the exclusive Australian distributor for companies including, but not limited to Mercury Medical, Allied Healthcare Products, VeinLite and Oxylitre.



NEANN MEDICAL KITS

Why buy from us? NEANN is Australia's largest manufacturer and leading designer of emergency medical kits supplying soft packs to 90% of Australia's Ambulance Services.

Our Medical Kits are manufactured from the most modern, state of the art material. Historically, our kits have been manufactured using Vironfrei material. Our kits have evolved through research and development to be produced using our new material called EnvironFrei. EnvironFrei has been independently tested, confirming its superiority to common place materials ensuring a stronger, longer lasting and more durable product.

At NEANN we understand the importance of cross-contamination and infectious control issues in the workforce. EnvironFrei directly combats this every increasing issue. Innovation in this area is constantly being achieved to deliver the best results to

meet our customers needs. A recent example of this would be our aVector handles which are designed to eliminate the effects of cross-contamination on our medical kits at the highest area of risk.

Our products are designed so that everything has a place and everything is in its place, there are no search and select issues with NEANN products. Layout guides and access panels are supplied with most NEANN products to aid the user with one of the most important issues, saving time.

Do you know what are you looking for in your kit? Do you require your kit to be tailored to your working needs? Do you know that NEANN can help you customise a kit or bag to any budget that will suit you or your business? With years of design experience, NEANN can deliver a product to suit your needs, in addition we can alter any of the bags from our range above.

RESPIRATORY & TRAUMA RANGE

air-Q sp Intubation Laryngeal Airway (ILA)

The Only Self-Pressurizing Mask Cuff.

Takes the guesswork out of mask cuff inflation because it:

- Eliminates the extra step of inflating the mask cuff,
- Eliminates mask cuff over inflation.

The new Self-Pressurizing air-Q sp enables intubation to be performed with standard size ET tubes through the device. Whilst the patient is effectively intubated you can remove the air-Q sp leaving the ET tube in place during the removal.



Product Code: 50-006-104035-0000



AG Cuffill

The AG Cuffill is designed to enhance patient safety by minimising the potential human error while measuring and inflating airway cuffs, as well as diminishing the risk of cross-contamination, when used for a single patient.

This pocket-size, syringe-like design allows a simple and easy operation and is used by medical professionals, including first-responders and hospital staff.

The compact and intuitive AG Cuffill device is the most accurate solution for measuring both pressure and volume of airway cuffs in all clinical settings.

Product Code: 50-063-HSCUFF0041-0000

air-Q sp / Blocker Intubation Laryngeal Airway (ILA)

Just Released—the NEW air-Q sp / Blocker Range.

Providing the ideal solution in one air-Q sp / Blocker airway that allows:

- Rescue ventilation,
- Rescue intubation,
- Rescue suctioning and venting of the esophagus.

The air-Q sp / Blocker is a new design incorporating a self-pressurising mask with a soft, flexible guide tube allowing access to the posterior pharynx and upper esophagus.



Product Code: 50-006-10-7045-0000



Product Code: 50-006-JQ-NF206-0000

naso-flo Nasopharyngeal Airway Device

The naso-flo Nasopharyngeal Airway Device allows for direct oxygen tubing connection and humidification vents for heat and moisture transfer on inhalation.

- Oxygen port and elbow allow for direct connection to O₂ source (higher FiO₂ % delivery),
- Soft silicone material designed for less traumatic insertions,
- Humidification vents allow for heat/moisture transfer on inhalation,
- Colour-coded elbows aid in correct size selection,
- Optional filtered respiratory indicator (CO₂ monitor connection).

RESPIRATORY & TRAUMA RANGE



RES-Q-VAC

RES-Q-VAC is a high-quality non-electric suction pump with removable single patient use canisters and catheters.

Reliable hand-operated advanced suction system which provides a full range of suction options for every medical requirement.

- Convenient one-handed operation,
- Portable and power safe,
- Always ready to go.

Product Code: 50-018-P1N1S000N-0000

CPR-2 Disposable Bag

Positive Pressure Ventilation and Self-Inflating Bags.

The CPR-2 range is a full line of resuscitators manufactured with new thermoplastic bag material providing superior lung compliance “feel” with exceptional responsiveness.

Benefits:

- Colour-coded Manometer, accuracy is $\pm 3\text{cm H}_2\text{O}$ between 0-15cm H_2O ,
- Monitors both Airway and PEEP Pressure,
- Fast re-expansion/high breath-per-minute capability,
- Latex-Free medical grade thermoplastic bag material.



Product Code: 50-006-1056002-0000



Full Range of CO₂ Detectors

Fast, durable colorimetric breath indicators for visualisation of exhaled CO₂ to assist in verifying proper intubation.

Features:

- Convenient Pull Tab activates the device when ready to use,
- Detects up to 50 breaths per minute,
- 24-hr performance,
- Works in 100% humidity,
- Vivid breath-to-breath colour change,
- Suitable patient weight starting from 0.25kg.

Product Code: 50-006-10-55370-0000

Manometer

The first disposable manometer to reliably monitor BOTH proper inflation pressure and PEEP pressure during manual ventilation.

Features

- The Disposable Colour-Coded Manometer allows for in-line monitoring of both airway pressure and PEEP pressures, when used in conjunction with CPR and hyperinflation resuscitation bags.
- The Disposable Colour-Coded Manometer is ready when you need it so you don't need to search or deal with a broken reusable manometer again!
- Disposable & Latex Free.



Product Code: 50-006-1055350-0000

RESPIRATORY & TRAUMA RANGE

Flow-Safe CPAP II - Disposable A Revolution in Emergency Care

- When time is of the essence,
- When accuracy is critical.

Features:

- Over 50% Less Oxygen Consumption. It uses less oxygen while delivering high FiO₂,
- Built-In Manometer. It verifies delivered CPAP pressure,
- Consistent CPAP Pressure,
- Very comfortable head harness and face mask.

**Implemented by
multiple ambulances
services in
Australia**



Product Code: 50-006-1057209-0000



Product Code: 50-092-SN-02L-0192

Ampoule Openers

Health Care Industries use millions of ampoules annually. Breakages, cuts and contamination of contents have always been a safety concern resulting in wasted resources.

The Snapit Solution is simple, safe and cost effective. Snapit is made from durable, anodised aluminium ensuring a once only investment provides ongoing safety.

- Protects the hands,
- An effortless open every time due to extended lever action,
- "Just snap it." It's a reusable device - hence great value!

VeinLite LEDX

VeinLite LEDX's C-Shaped design and 32 LEDs in two colours means that superficial veins can be seen and accessed anywhere on the body. It secures the vein in place and stretches the skin for easy needle insertion without any guess work.

Features:

- Best ultra portable system on the market,
- Disposable plastic covers protect the device from contamination & cross-contamination,
- Light shield and paediatric adapter supplied,
- Weight including battery - 83g.



Product Code: 50-019-VLEDX-0000



Product Code: 50-001-0734-0000

NEANN LSB (Long Spine Board)

This non-conductive, super light NEANN Long Spineboard (LSB) delivers superior comfort for paramedics with newly designed hand holds, curvature angles and increased strength, making manual handling easier while still maintaining the functionality of the spineboard.

Features:

- Light weight,
- Greater flexibility,
- Safety tested.

MEDICAL & FIRST AID KIT RANGE

The NEANN Professional Range of Kits are the latest NEANN designs, utilising the state of the art manufacturing methods and the newly developed EnvironFrei materials and aVector handles, this is the kit range of the year!

NEANN kits achieve their pre-eminence in the Australia market with meticulous attention to detail and an unrelenting pursuit of continual improvement.



Product Code: 50-001-0171-0192

Hospital In The Home (HITH) Equipment Kits

Specifically designed for the Hospital in the Home (HITH) professional the Equipment Kit will carry the OBS equipment and consumables required for a clinician to provide care for 8 to 10 daily visits. The accompanying HITH Vehicle Kit is designed to carry a laptop and patient notes as well as patient specific drugs and fluids. The HITH Vehicle Kit is easily attached to a folding trolley for ease of transport.

Features:

- Light and convenient with replaceable parts,
- Product specific allocation including personal items.

Drug IV Professional Series Kits

The Drug IV Professional Series Kits have been design with equipment laid out in order of use. This assists the user in establishing IV access with ease. This Drug IV Professional Series Kit is design to fit into the PRK Professional Series Kit.

Features:

- Stores large range and number of ALS drugs,
- Unique drug labelling system for easy identification of drugs.



Product Code: 50-001-8300-0386

OxyRESQ

With 2 compartments - Top Pouch and Main Oxygen Therapy Bag the OxyRESQ is designed to be flexiabe and carry the equipment you want, to the level of care you wish to provide. A full range of basic to advanced airway equipment can be housed in this kit.

Features:

- Lightweight and ergonomically designed,
- Houses our new disposable venturi suction unit.



Product Code: 50-001-1063-0369

TEK Professional Series Kits

The TEK Professional Series Kits provides a place for everything and everything in its place with removable wear base, heavy duty runners and highly reflective tape for safety, this is a highly durable kit with great features.

Features:

- Easy to read printed labelling,
- Built-in layout guides.



Product Code: 50-001-1455-0192

WHY PURCHASE OUR NEANN KITS?

Here is an overview of why NEANN kits continue to lead the market from our competitors in our pursuit of producing Effective, Reliable and Safe to Use products:

We can Design, Research & Develop a kit / bag for you!

Design the kit to the level of care and equipment you want to carry and suit the environment, emergency or everyday situations you may come across while in the workforce.

EnvironFrei

Modern State of the art material.
Stronger.
Highly Durable.
Improved Infectious Control.
Sustainability.

Cross-Contamination

Unique material design to minimise Cross-Contamination of bacteria and improved infectious control issues in the workforce.
Easy to clean off many different types of bacteria, infections and grime.

Product Labelling

No Search and Select Issues.
Unique internal labelling system for rapid location and easy identification of products in an emergency.
Easy to re-stock.

Handles & Straps

aVector handles.
Improved Infectious Control.

Clear Open Design

Allows all gear to be easily located.
Everything has a place and everything in its place.

Independent Reporting & Testing

Performance Based.
Ergonomically Tested & Approved.

Layout Guides

To assist user in quickly identifying and re-stocking products.
Saving Time.

Logo Design and Placement

Your unique logo and identification (for example: company name, address) can be added to the kit.



Modular System

So you can slowly build up the kits as finances become available.
Internal and External Fixed Storage Pockets also available.

Removable Panels / Inserts

Internal removable lid panels can be interchangeable to suit multiple uses when on the run or if situations change.

Runners & Base

Heavy duty runners and removable wear base thus improving the life of the product.

Soft Packs

Padded with stiffener and foam, these kits are design to endure the harsh outdoor environments while being soft and padded for easy storage.

Multiple Colour Options For Your Kit

Colour co-ordinate the material for the type of equipment housed.
Light and easy to see colours.

Multiple External Window Panels

For instructions, manual and layout guides to assist users and for increased safety.

Reflective Tapes & Zips

High visibility & safety.

Padded Backpack

Designed ergonomically to provide comfort and support in all terrains.

Separate Dividers or Elastic Holders

To keep products in their place.
Prevent spillages.
Ensuring equipment is always secure.

TGA Certified

Our kits are Certified to ensure to you are getting the best quality product that complies with the Australian Therapeutic Goods Administration (TGA).

Warranty

All kits supplied should be free from defect, should this not be the case we will replace or repair the item (at our option).



For more information please contact
RAPP Australia Pty Ltd

PO Box 120
Lara, Victoria 3212, Australia

Office: +61 3 5284 0222
Fax: +61 3 5282 4892

Email: sales@rappaustralia.com.au
www.neann.com.au

